
AX224 Biological Indicator

Rev.01.2025

Quality Certificate / Kalite Sertifikası

Lot	:	_____
<i>Geobacillus stearothermophilus</i>		
Manufacturing Date	:	_____
Üretim Tarihi	:	_____
Expiration Date	:	_____
Son Kullanma	:	_____
Population	:	_____ CFU/KOB
Popülasyon	:	_____
D-Value (121°C)	:	_____ min/dk
D-Değeri	:	_____
Survival Time	:	_____ min/dk
Hayatta Kalma Süresi	:	_____
Survival time = not less than (log labeled population - 2) x labeled D-value		
Kill Time	:	_____ min/dk
Ölüm Zamanı	:	_____
Kill time = not more than (log labeled population + 4) x labeled D-value		
D-Value (132°C)	:	_____ sec/sn
D-Değeri	:	_____
Survival Time	:	_____ min/dk
Hayatta Kalma Süresi	:	_____
Survival time = not less than (log labeled population - 2) x labeled D-value		
Kill Time	:	_____ min/dk
Ölüm Zamanı	:	_____
Kill time = not more than (log labeled population + 4) x labeled D-value		
D-Value (134°C)	:	_____ sec/sn
D-Değeri	:	_____
Survival Time	:	_____ sec/sn
Hayatta Kalma Süresi	:	_____
Survival time = not less than (log labeled population - 2) x labeled D-value		
Kill Time	:	_____ min/dk
Ölüm Zamanı	:	_____
Kill time = not more than (log labeled population + 4) x labeled D-value		
Z-Value	:	_____ °C
Z-Değeri	:	_____

AXIS IS A BRAND OF KORDON TIP LTD.

Manufacturer / Üretici

Kordon Tıp Sağlık Araç ve Gereçleri
Müh.Proje İth.San. ve Tic.Ltd.Şti.
10006/1 No:43 AOSB Çiğli / İzmir TÜRKİYE
Tel: +90 232 348 60 60
Fax: +90 232 348 60 90
www.kordontip.com - info@kordontip.com

EN AX224 Steam Ultra Rapid Biological Indicator

Product Description

AXIS® AX224 Steam Ultra Rapid Biological Indicator is specially designed for monitoring steam sterilization processes. If the sterilization process is unsuccessful, after incubation at 60±2°C, the medium will give a positive fluorescent readout from the **AXIS®** Automatic Reader Incubator, indicating the presence of viable *Geobacillus stearothermophilus* spores. After incubation, if the sterilization process is successful, the **AXIS®** Automatic Reader Incubator will not show a fluorescence readout and will have a negative result. The last readout should be taken 20 minutes after incubation at 60±2°C.

Content

AXIS® AX224 Steam Ultra Rapid Biological Indicator consists of medium in a glass ampoule, a carrier inoculated with *Geobacillus stearothermophilus* spores, a chemical indicator strip, a plastic tube, a plastic cap, and a filter paper inside the cap.

Incubation of **AXIS®** AX224 Steam Ultra Rapid Biological Indicators should be done in the **AXIS®** Automatic Reader Incubator. The fluorescent readout is produced when the reader stimulates the medium with UV light at 360 nm. The fluorescent readout is an indirect detection of the germination and growth of *Geobacillus stearothermophilus* spores that survive in **sufficient population** after sterilization. A failure in the sterilization process may also be detected by a color change in the medium. Incubation of **AXIS®** AX224 Steam Ultra Rapid Biological Indicators longer than 20 minutes is not an advantage due to the high sensitivity of the fluorescent readout.

Instructions for Use

1. **AXIS®** AX224 Steam Ultra Rapid Biological Indicators have an identification label along the tube where the user can write details such as sterilizer number, load number and date of operation.
 2. Place the biological indicators together with the materials to be sterilized in separate packages according to the recommended sterilization practices. Place this package in areas where you think the sterilizing agent (steam) will have the hardest reach. These difficult areas are generally the center of the load and areas near the door.
 3. Perform your routine sterilization process.
 4. After the sterilization and aeration process is finished, wear protective glasses and gloves and take out the **AXIS®** AX224 Steam Ultra Rapid Biological Indicator from the sterilized package, and wait 5-10 minutes until it cools down to room temperature.
 5. Check the chemical indicator strip on the biological indicator label. The absence of any color change in the chemical indicator indicates that the biological indicator has not been exposed to the sterilizing agent. A color change from pink to brown or its shades indicates exposure to the sterilizing agent.
- Note:** The color change in the chemical indicator does not guarantee sterilization.
6. Break the glass ampoule in the biological indicator using the ampoule breaking socket located in the incubation area of the incubator. After breaking the glass ampoule, make sure that the spore paper is completely wetted by the medium, and immediately place the **AXIS®** AX224 Steam Ultra Rapid Biological Indicator into the **AXIS®** Automatic Reader Incubator.
7. If the sterilization process is unsuccessful, the medium will still give a positive fluorescent readout from the **AXIS®** Automatic Reader Incubator after incubation at 60±2°C, indicating the presence of still viable *Geobacillus stearothermophilus* spores. If the sterilization process is successful, after incubation, the **AXIS®** Automatic Reader Incubator will not read fluorescence and will have a negative result. The last readout should be taken 20 minutes after incubation at 60±2°C. If there is no color change in the medium of the indicator exposed to the sterilization process at the end of 20 minutes, and if the medium has not been exposed to the sterilization process and the color of the medium of the positive control indicator has turned yellow or shades, the sterilization process has been completed successfully. Record the results and dispose of any positive biological indicators immediately as described below.

Note: When performing routine sterilization, a nonsterile biological indicator should be used as a positive control at least once a day. A positive control indicates correctness of incubation conditions, spore viability, adequacy of the medium, and that the **AXIS®** Automatic Reader Incubator is reading correct fluorescence. For the positive control indicator result to be valid, a fluorescent reading must be taken from the **AXIS®** Automatic Reader Incubator within 20 minutes at 60±2°C. The positive control indicator and the processed indicator must belong to the same lot.

Visual Confirmation: 24 hours

Optionally, the success of the sterilization process can be achieved by visual color change tracking after 24 hours of incubation. If the sterilization process is unsuccessful, after incubation at 60±2°C, the medium will change from purple to yellow or its shades, indicating the presence of still viable *Geobacillus stearothermophilus* spores. If the sterilization process is successful, there will be no change in the color of the medium after incubation, it will remain in its original color. If incubation is longer than 24 hours, a moist environment is required to prevent the medium from drying out.

Storage

Keep the products in their original boxes and store them at a temperature of 10-30°C and a relative humidity of 30-80%. Do not freeze. Indicators should not be stored near sterilants or other chemicals. **AXIS®** AX224 has a shelf life of 24 months.

Disposal

Dispose of the biological indicator in accordance with the health and safety procedures applicable in your country. A positive biological indicator can be autoclaved at 121°C for 30 minutes, at 132°C for 15 minutes or at 134°C for 10 minutes in a gravity steam sterilizer; and at 132°C for 4 minutes or at 135°C for 3 minutes in a high speed pre-vacuum sterilizer.

Warnings

I Do not use **AXIS®** AX224 Steam Ultra Rapid Biological Indicator for EtO, dry heat, formaldehyde, hydrogen peroxide and other sterilization processes.
! Do not reuse the biological indicators.

AXIS® AX224 Steam Ultra Rapid Biological Indicator is in compliance with **ISO 11138-1:2017** and **ISO 11138-3:2017**.

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