
# AX100 Biological Indicator

Rev. 08.2023

## Quality Certificate / Kalite Sertifikası

Lot :  
Geobacillus stearothermophilus (ATCC 7953)

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Expiration Date :  
Son Kullanma

Population : CFU/KOB  
Popülasyon

D-Value : min/dk  
D-Değeri  
60°C, FORM 1mol/L

Survival Time : min/dk  
Hayatta Kalma Süresi  
Survival time = not less than (log labeled population - 2) x labeled D-value

Kill Time : min/dk  
Ölüm Zamanı  
Kill time = not more than (log labeled population + 4) x labeled D-value

AXIS IS A BRAND OF KORDON TIP LTD.

### Manufacturer / Üretici

Kordon Tıp Sağlık Araç ve Gereçleri Ltd. Şti.  
10006/1 No:43 AOSB Çiğli / Izmir TÜRKİYE  
Tel: +90 232 348 60 60  
Fax: +90 232 348 60 90  
www.kordontip.com - info@kordontip.com

# EN AX100 Formaldehyde Biological Indicator

## Product Description

**AXIS®** AX100 Formaldehyde Biological Indicator is specially designed for monitoring ethylene oxide sterilization processes. If the sterilization process is unsuccessful, after incubation at 60±2°C, the color of the medium will change from green to yellow or its shades, indicating the presence of still viable *Geobacillus stearothermophilus* (ATCC 7953) spores. If the sterilization process is successful, there will be no change in the color of the medium after incubation, it will remain in its original color. The last readout should be taken after 48 hours of incubation at 60±2°C.

## Content

**AXIS®** AX100 Formaldehyde Biological Indicator consists of medium in a glass ampoule, a disc inoculated with *Geobacillus stearothermophilus* (ATCC 7953) spores, a chemical indicator strip, a plastic tube, a plastic cap, and a filter paper inside the cap.

## Instructions for Use

1. **AXIS®** AX100 Formaldehyde Biological Indicators have an identification label along the tube where the user can write details such as sterilizer number, load number and date of operation.
2. Place the biological indicators together with the materials to be sterilized in separate packages according to the recommended sterilization practices. Place this package in areas where you think the sterilizing agent will have the hardest reach. These difficult areas are generally the center of the load and areas near the door.
3. Perform your routine sterilization process.
4. After the sterilization and aeration process is finished, wear protective glasses and gloves and take out the **AXIS®** AX100 Formaldehyde Biological Indicator from the sterilized package, and wait 5-10 minutes until it cools down to room temperature.
5. Check the chemical indicator strip on the biological indicator identification label. The absence of any color change in the chemical indicator indicates that the biological indicator has not been exposed to the sterilizing agent. A color change from pink to brown or its shades indicates exposure to the sterilizing agent.
- Note:** The color change in the chemical indicator does not guarantee sterilization.
6. Break the glass ampoule in the biological indicator using the ampoule breaking socket located in the incubation area of the incubator. After breaking the glass ampoule, make sure that the spore paper is completely wetted by the medium, and incubate the **AXIS®** AX100 Formaldehyde Biological Indicator at 60±2°C as soon as possible.
7. If the sterilization process is unsuccessful, after incubation at 60±2°C, the medium will change from green to yellow or its shades, indicating the presence of still viable *Geobacillus stearothermophilus* (ATCC 7953) spores. If the sterilization process is successful, there will be no change in the color of the medium after incubation, it will remain in its original color. The last readout should be taken 48 hours after incubation at 60±2°C. If there is no color change in the medium of the indicator exposed to the sterilization process at the end of 48 hours, and if the medium has not been exposed to the sterilization process and the color of the medium of the positive control indicator has turned yellow or shades, the sterilization process has been completed successfully. Record the results and dispose of any positive biological indicators immediately as described below.
- Note:** An unsterilized biological indicator should be used as a positive control at least once a day when performing routine sterilization. A positive control indicates the correctness of the incubation conditions, spore viability, and adequacy of the medium. For the positive control indicator to be valid, it should turn yellow or shades of yellow after 48 hours of incubation at 60±2°C. The positive control indicator and the processed indicator must belong to the same lot.

## Storage

Keep the products in their original boxes and store them at a temperature of 10-30°C and a relative humidity of 30-80%. Do not freeze. Indicators should not be stored near sterilants or other chemicals. **AXIS®** Ax100 has a shelf life of 24 months.

## Disposal

Dispose of the biological indicator in accordance with the health and safety procedures applicable in your country. A positive biological indicator can be autoclaved at 121°C for 30 minutes, at 132°C for 15 minutes or at 134°C for 10 minutes in a gravity steam sterilizer; and at 132°C for 4 minutes or at 135°C for 3 minutes in a high speed pre-vacuum sterilizer.

## Warnings

! Do not use **AXIS®** AX100 Formaldehyde Biological Indicator for steam, dry heat, hydrogen peroxide and other sterilization processes.  
! Do not reuse the biological indicators.

**AXIS®** AX100 Formaldehyde Biological Indicator is in compliance with **ISO 11138-1:2017** and **ISO 11138-2:2017**.

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