

# BIOLOGICAL INDICATOR COLOR CONVERSION CHART



**AX224 - 20 min Steam BI**

## CHEMICAL INDICATOR

### *PRIOR STERILIZATION*

**AXIS**

AX 224

*Geobacillus stearothermophilus*

Ster. | Load | Date

**Lot**

### *AFTER STERILIZATION*

**AXIS**

AX 224

*Geobacillus stearothermophilus*

Ster. | Load | Date

**Lot**

## CULTURE MEDIUM

*Negative(-)*



*Positive(+)*



THE COLOR ON THE PRODUCT CAN SHOW DIFFERENCE WITH THE COLOR SEEN ON REFERENCE CHART DUE TO PRINTING QUALITY. USER CAN MAKE COMPARISON WITH UNUSED PRODUCT.