



Chemical Indicators Solutions & Features

Product Code	Indicator Type	Cycle	Initial Color Final Color	Adhesive
AXCP90	Type 5	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✗
AXCP91	Type 5	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✗
AXCP350	Type 5	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✓
AXCP60	Type 6	3.5 min @ 134 °C	 → 	✗
AXCP61	Type 6	5 min @ 134 °C	 → 	✗
AXCP70	Type 6	7 min @ 134 °C	 → 	✗
AXCP450	Type 6	3.5 min @ 134 °C	 → 	✓
AXSDC3	Type 4	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✗
AXSDC4	Type 4	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✓
AXSDC5	Type 4	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✗
AXSDC6	Type 6	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✓
AXSDC7	Type 4	3.5 min @ 134 °C / 15 min @ 121 °C	 → 	✗