

@x **AXIS**[®]

TYPE 2



compatible with ISO 11140-1:2017
ISO 11140-1:2017 ile uyumludur.



Kordon Tıp Sağlık Arac ve Gereçleri Ltd.
10006/1 No:43 AOSB Cigli IZMIR TURKEY
Tel:+90 232 348 60 60
Fax: +90 232 348 60 90
www.kordontip.com - info@kordontip.com

Discard after use, according to your country's healthcare and safety regulations.

Kullandıktan sonra, ülkenizin sağlık ve güvenlik düzenlemelerine göre yok edin.

Figure 1 illustrates the color calibration process. It shows four rows of color patches, each labeled with 'PCD-4S', 'STEAM', and '60x AXIS'. The rows are categorized as follows:

- Unprocessed:** The first row shows a series of pink squares.
- Failures- any color rather than reference:** The second row shows a series of brown and pink squares.
- Processed and satisfactory:** The third row shows a series of brown squares.
- Processed and satisfactory:** The fourth row shows a series of black squares.

ENG

	PCD-4S	STEAM	AXIS	İşlem görmemiş
	PCD-4S	STEAM	AXIS	Hatalı Sonuç - referans rengi hariç tüm renkler
	PCD-4S	STEAM	AXIS	İşlem görmüş ve başarılı
	PCD-4S	STEAM	AXIS	

TR